

\*  $\mu\text{L} \rightarrow \text{uL}$ !

more than samples  $\rightarrow$  96 well plate

\* can do sample prep night before + put homogenates in  $-80^\circ\text{C}$ , till next day.

## RNA Extraction with RNA SPRI Beads:

\* make sure you have enough of each buffer

sample prep

24 or thicker area facing outwards

\* take a photo of order of tubes

if marker may rub off in tissue lyzer

\* make it!

mark of tubes on each side

Take pic of tubes in rack

so you know

what sample

is in each well.

check

1. Pull out and place in a tube the amount of RNA beads you need ( $90\mu\text{L}/\text{sample}$ ;  $n+1$  amount) to allow them to acclimate to room temp.  $90 \cdot n + 10\%$  (take  $10\%$  of  $n$  and add it)
2. Place each sample in a 2mL tube with 1 clean BB.
3. Add 100uL lysis buffer to each sample. Shake the tubes in the Tissue Lyzer for 1-3min at 30hz.
4. Add 60uL 100% isopropanol to each tube. Mix by pipetting and incubate for 1 minute.
5. Spin samples down in the centrifuge for 1 min at  $12,000g \rightarrow \text{RCF NOT RPM}$
6. In a 96-well PCR plate (many samples) OR in autoclaved 1.5mL tubes (few samples), add 90uL RNA magnetic beads to each well/tube and 10uL Lysis/Binding Enhancer (LBE). Mix well.  $\rightarrow$  centrifuge
7. Move the samples from centrifuge tubes to either the plate or tube. Pipette up and down 15-20 times to mix the sample and beads together well. Incubate at room temp for 5-10 minutes. (binds to beads)  $\rightarrow$  try not to get floaties.
8. Place the plate/tubes on a magnetic rack and allow the beads to separate (2-3min). Remove the supernatant. NOTE: The buffers in this step cause the beads to come off the magnet easily, so be gentle when removing supernatant to not remove beads.  $\rightarrow$  take off rack
9. Pull plate/tubes off the magnetic rack. Add 150uL RNA Beads Wash Buffer (WB) 1, pipette up and down 10-15 times to mix the buffer and beads well. Incubate at room temp for 2 min.  $\rightarrow$  no bubbles!  $\rightarrow$  yellow tip
10. Place the plate/tubes on the magnetic rack and allow the beads to separate (2-3min). Remove the supernatant being careful not to remove beads.  $\rightarrow$  use multichannel set @  $150\mu\text{L}$
11. Pull plate/tubes off the magnetic rack. Add 150uL WB2, pipette up and down 15-20 times to mix the buffer and beads well. Incubate at room temp for 2 min.  $(n+10\%) \rightarrow$  reservoir  $(150 \cdot n + 10\%)$
12. Place the plate/tubes on the magnetic rack and allow the beads to separate (2-3min). Remove the supernatant being careful not to remove beads. During the separation of beads and WB2, make a mixture of DNase and DNase buffer for your samples. For each sample add 24uL WB2, 3uL DNase and 3uL DNase Buffer. Make  $n+1$  amount.  $\rightarrow n+10\%$  amount  $\rightarrow$  separate?
13. After removing WB2 supernatant, remove the plate/tubes from the magnet. Add 30uL of the DNase mixture to each sample. Pipette up and down 10-15 times to mix beads and DNase together. Incubate at room temperature for 15-30min.
14. After incubating, add 100uL of the binding buffer to each of the samples. Pipette up and down to mix well. Incubate for 3-5min.  $(n+10\%)$
15. Place the plate/tubes on the magnetic rack and allow the beads to separate (2-3min). Remove the supernatant being careful not to remove beads.  $\rightarrow n+10\%$
16. Pull plate/tubes off the magnetic rack. Add 150uL WB2, pipette up and down 15-20 times to mix the buffer and beads well. Incubate at room temp for 2 min. Place the plate/tubes on the magnetic rack and allow the beads to separate (2-3min). Remove the supernatant being careful not to remove beads. Repeat this for a total of 2 times. (steps 15 + 16)
17. On the second removal of WB2, leave the plate/tube on the magnetic rack. Let the beads air dry for 3-5min, but being careful not to over dry them (over drying causes lower yield).  $\rightarrow$  closer to 3  $\rightarrow$  on magnetic rack
18. Remove the plate/tubes from the magnetic rack after drying. Add 25-50uL autoclave ddH<sub>2</sub>O to each sample to elute the RNA. Pipette up and down 10-15 times to mix the beads and water well. Incubate at room temp for 5 min.
19. Place the plate/tubes on the magnetic rack and allow the beads to separate (2-3min). Remove all but 2uL of the water from the plate/tubes being careful not to remove beads. Place the RNA in new tubes.
20. Qubit the RNA for concentration and store at  $-80^\circ\text{C}$  for downstream experiments. Take off 23uL

- Any time something goes into hood, spray w/ ethanol.
- spray gloves down if you leave hood.
- can set timers on phone - gloves off, use, new gloves

\*  $n+10\%$  only in steps.

$100\mu\text{L} \times n + 10\% = \text{needed volume}$

MOSSQUITO